

Dyename HT-holder for PEC-measurements

- Short Description -

Small holes for degassing or gas extraction for analysis

Screws can be modified to combine with other instruments


Diameter: 2 mm

Hole for additional electrode (e.g. reference electrode)

Can be sealed by septum or solid screw

Can be modified to combine with other instruments

Diameter: 6 mm



Place for membranes

Nafion 211 membrane is included

Diameter of window: 12 mm

Space for photoelectrode

(between Teflon and metal mask)

Diameter of light window: 8 mm